Modern
Med

EUROCOLUMBUS ALIEN X

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EUROCOLUMBUS ALIEN X

Aparat RTG z ramieniem C

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Gwarancja do 24.10.2024

Moc generatora 25kW

Detektor 30x30 cm

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ZBIGNIEW ROMANOWSKI

+48 695 543 800

zbignew.karo@gmail.com

ALIEN-X

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A) Unit general conditions

N.	Before maintenance	Maintenance check/procedure	After maintenance
1	☒ PASS ☐ FAIL	Check integrity of the external covers	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check integrity of the external cabling (e.g. displays, accessories)	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check integrity of the umbilical connector, both mechanical and electrical and optical	☒ PASS ☐ FAIL
4	☒ PASS ☐ FAIL	External cleaning check	☒ PASS ☐ FAIL
5	☒ PASS ☐ FAIL	Internal cleaning check PASS	☒ PASS ☐ FAIL

B) Mechanical functions

N.	Before maintenance	Maintenance check/procedure	After maintenance
1	☒ PASS ☐ FAIL	Check integrity of the controlling handles	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check steering and brakes	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check manual movement devices (wig-wag, horizontal carriage)	☒ PASS ☐ FAIL
4	☒ PASS ☐ FAIL	Check trolley wheels and their brakes	☒ PASS ☐ FAIL
5	☒ PASS ☐ FAIL	Check the C-arm manual motion activation devices (1)	☒ PASS ☐ FAIL

(1) only on unit equipped with manual motorized C-arm

C) Electromechanical functions

N.	Before maintenance	Maintenance check/procedure	After maintenance
1	☒ PASS ☐ FAIL	Check orbital motion and its position sensor	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check circumferential motion and its position sensors and end sensors	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check horizontal carriage brake, check the horizontal motor if present and horizontal position sensor if present (2)	☒ PASS ☐ FAIL
4	☒ PASS ☐ FAIL	Check lift column and its position sensor	☒ PASS ☐ FAIL
5	☒ PASS ☐ FAIL	Check DFF (detector lift) and its position sensor	☒ PASS ☐ FAIL
6	☒ PASS ☐ FAIL	Check trolley monitors lift actuator	☒ PASS ☐ FAIL
7	☒ PASS ☐ FAIL	Check motion of the trolley touch control panel	☒ PASS ☐ FAIL

(2) only on unit equipped with motorized horizontal carriage

(3) only on unit equipped with manual/motorized C-arm

D) X-ray inverter and generator tank

N.	Before maintenance	Maintenance check/procedure	After maintenance
1	☒ PASS ☐ FAIL	Check HT connector, integrity, clean and refill silicon grease.	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check emission calibration and even recalibrate filaments.	☒ PASS ☐ FAIL

3	☒ PASS ☐ FAIL	Check dose linearity, constancy. (*)	☒ PASS ☐ FAIL
---	--	--------------------------------------	--

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E) tubehead and collimator

N	Before maintenance	Maintenance check/procedure	After maintenance
.			

1	☒ PASS ☐ FAIL	Check HT connector, integrity, clean and refill silicon grease	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check liquid cooling system	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check tube efficiency. (*)	☒ PASS ☐ FAIL
4	☒ PASS ☐ FAIL	Check collimator positioning and even recalibrate it.	☒ PASS ☐ FAIL

(*) usually performed by local radioprotection service

F) Panel detector

N	Before maintenance	Maintenance check/procedure	After maintenance
.			

1	☒ PASS ☐ FAIL	Detector recalibration.	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check laser centering tool and even recalibrate	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check field illumination lamp, if present (5)	☒ PASS ☐ FAIL

(5) optional

G) Software

N	Before maintenance	Maintenance check/procedure	After maintenance
.			

1	☒ PASS ☐ FAIL	Check turning on time	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check graphic user interface functionality	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check dose point reproducibility (AEC)	☒ PASS ☐ FAIL
4	☒ PASS ☐ FAIL	Image quality test	☒ PASS ☐ FAIL
5	☒ PASS ☐ FAIL	X beam alignment and uniformity test	☒ PASS ☐ FAIL
6	☒ PASS ☐ FAIL	Check functioning of wireless pedal and remote controls, if present (6)	☒ PASS ☐ FAIL

(6) optional

H) Monitors

N	Before maintenance	Maintenance check/procedure	After maintenance
.			

1	☒ PASS ☐ FAIL	Check integrity of the image displays	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check functioning and even calibrations of the image displays	☒ PASS ☐ FAIL

S) Sicurezza

N	Before maintenance	Maintenance check/procedure	After maintenance
.			

1	☒ PASS ☐ FAIL	Check functioning of the emergency switch.	☒ PASS ☐ FAIL
2	☒ PASS ☐ FAIL	Check functionality of the xray interlock key switch.	☒ PASS ☐ FAIL
3	☒ PASS ☐ FAIL	Check integrity and functionality of the xray emission alarms, both visual and acoustic.	☒ PASS ☐ FAIL
4	☒ PASS ☐ FAIL	Check integrity and functioning of the service port (injector, doorswitch, external exposure lamp) if used/present (7)	☒ PASS ☐ FAIL
5	☒ PASS ☐ FAIL	Check functioning of the main protection breaker.	☒ PASS ☐ FAIL

6	☒ PASS ☐ FAIL	Check/inspection the main power plug, especially the ground connection	☒ PASS ☐ FAIL
7	☒ PASS ☐ FAIL	Perform the electrical safety test according with the IEC 60601-1 and collateral directives.	☒ PASS ☐ FAIL

(7) optional

Notes about the performed corrections:

(please indicate the maintenance step, for example H7 if related to the final electrical safety test, and provide details about the activity performed in order to pass the maintenance point)

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