

# PDF

-25°C -40°C

## SUPER DEEP FREEZERS



PDF 270

PDF 370 - 440 - 530 - 625

PDF 370 - 440 - 530 - 625

MODEL	270	370	440	530	625
Capacity (litres)	270	370	440	530	625
Dimensions (W x D x H cm)	60 x 76 x 167	60 x 64 x 192	60 x 74 x 192	75 x 74 x 210	75 x 86 x 210
External structure	W / S	W / S	W / S	W / S	W / S
Internal chamber	W / S	W / S	W / S	W / S	W / S
Insulation	50	50	50	75	75
Feet		4	4	4	4
Rollers	4				
Insulated door	1	1	1	1	1
Internal lighting	•	•	•	•	•
Shelf	2	3	3	3	3
Digital electronic thermostat with alarm	•	•	•	•	•
Control panel xPRO	•	•	•	•	•
Forced-air refrigeration	•	•	•	•	•
Automatic defrosting	•	•	•	•	•
Automatic evaporation of the condensate water	•	•	•	•	•
Refrigerant	R404A	R404A	R404A	R404A	R404A
CFC- and HCFC- free	•	•	•	•	•
Voltage	230/1/50	230/1/50	230/1/50	230/1/50	230/1/50

W=White S=Stainless steel

ACCESSORIES	270	370	440	530	625
Shelf	•	•	•	•	•
Stainless steel wire shelf				•	•
Wire basket				•	•
Drawer	•	•	•	•	•
Drawer dividers	•	•	•	•	•
Drawer front top	•	•	•	•	•
Castors	•	•	•	•	•
Door key locking	•	•	•	•	•
Temperature chart recorder	•	•	•	•	•
Electronic temperature printer	•	•	•	•	•
Data-logger function with USB port	•	•	•	•	•
Ballasted sensor for simulation of product	•	•	•	•	•
GSM/GPRS telephone dialer	•	•	•	•	•
Voltage 240V/1ph/60Hz	•	•	•	•	•