

Silent PRIME versions with air dryer









CS 50/25 PRIME M








Code: PR-CS300-040-M

Description: Silent oil free compressor, FAD at 5 bar 225 l/min, tank capacity 40 lt, with membrane air dryer able to warranty till -20° C dew point value.

This compressor is available in 10 BAR VERSION also! visit the "10 bar oil-free compressors section" for more details!

Model	 Dental chairs	 Air dryer	 Air Capacity	 FAD @ 5bar	 Motor	 Volt.	 Freq.	 Receiver
			m ³ /h l/min cfm	m ³ /h l/min cfm	HP KW	V	Hz	l GAL
CS 50/25 PRIME M	4	yes	- 360 13,5	- 225 7,9	3 2,2	230	50	40 10.57

Model	 Noise Level	 Duty cycle	 Max pressure	 Size	 Net weight
	dB(A)		bar PSI	mm	Kg lb
CS 50/25 PRIME M	53	70%	7 101	490x720x890	106 233

