


EndoNaut®

- PAD -

Navigation, safety
and optimization.



- ✓ Less X-ray and contrast doses
- ✓ Optimized stent size
- ✓ Secured device deployment with augmented roadmaps



EndoNaut®

Designed to upgrade your OR at a lower cost

1 Easy generation of angiography panorama

- Clear and comprehensive view of lower limb






No OR / IT / HR constraints

Just plug EndoNaut® to the C-arm video output.
It's ready for use immediately!

2 Delimitation of lesions

- Lesion cues and length measurement



3 Real-time panorama-based roadmapping

- Lowers the need of contrast agent and X-rays
- Saves time during operations
- Compensates for table and patient movements
- Optimizes balloon / stent length sizes for better outcomes

🌐 www.Therenva.com
✉ contact@therenva.com

🏠 **Headquarter:**
74F rue de Paris
35000 Rennes, FRANCE
☎ +33 972 52 29 20

🏠 **US Office:**
761 Broadoak Loop
Sanford, FL 32771 USA
☎ +1 (407) 687-7054

EndoNaut® is a registered trademark of Therenva SAS.

In the USA, the EndoNaut® software obtained a substantial equivalence determination and FDA clearance through the CDRH premarket notification process (510(K)).

In Europe, the EndoNaut® software is CE marked (class IIb), not eligible for reimbursement.

The information provided in the labelling and manual is intended for healthcare professionals only. For the safe and successful operation and use of the device, always read the instructions. Update 2020/11/12. Non contractual document - ©2020 Therenva SAS. All rights reserved.