



# Confidence in your patient's comfort

Respironics patient interfaces for noninvasive ventilation

**PHILIPS**  
sense and simplicity



# Delivering optimal NIV

Respironics masks are designed to support enhanced NIV. From comfort for the patient to ease of use for the clinicians, our interfaces are specially designed to meet the rigors of acute care environments.



## Respironics PerformMax

### Key benefits

- Available for adult and pediatric patients (7 years or older and > 20kg)
- Perimeter facial contact increases comfort and simplifies fitting
- The clear, unobstructed view provided by the PerformMax total face mask helps reduce patient claustrophobia

Feature level



Compatible with a leak 2 entrainment elbow (EE), standard elbow (SE), and bronchoscopy elbow (BE)

## The first choice for comfort and fit

Total face masks allow quick fitting, and eliminate nasal bridge challenges by sealing around the perimeter of the face where patients have less pressure sensitivity and smoother facial contours.

### Mask selection criteria

- Rapid application
- Mouth breather
- Claustrophobia
- Facial abnormalities
- Lack of teeth

## Respironics Total

### Key benefits

- The clear, unobstructed view of the Respironics Total mask helps reduce patient claustrophobia
- Perimeter facial contact increases comfort and simplifies fitting
- Digital Auto-Trak compatibility enhances NIV performance

Feature level





## Respironics AF531

### Key benefits

- CapStrap headgear provides an excellent fit and simple re-application
- Support your infection control efforts with the CleanClip system
- You can adapt one mask to multiple ventilators with interchangeable elbows

Feature level



Compatible with a leak 1 or leak 2 entrainment elbow (EE), EE NIVO elbow (EE/NE), standard elbow (SE), SE NIVO elbow (SE/NE), and bronchoscopy elbow (BE)

## Taking comfort further

Oro-nasal masks are the standard of care for hospital NIV, and the Respironics family includes a comprehensive selection of systems with features that address patient comfort, mask fitting, and ventilator compatibility.

### Mask selection criteria

- Mouth breather
- Long-term NIV

## Respironics AF811

### Key benefits

- Cushion uses a gel pad with a silicone flap for increased comfort
- CapStrap headgear system simplifies mask application and routine patient care
- Available in three sizes to fit a wide patient population

Feature level



Entrainment elbow (EE)  
Also available with standard elbow (SE)



## Patient comfort made easy

CapStrap headgear provides excellent stability and simplifies initial setup by keeping the mask in place while you adjust the patient's straps. When you need to interrupt treatment for oral care and medications, the CapStrap headgear allows quick removal and simple re-application. CapStrap headgear is available with several Respironics masks:

- AF531
- AF811
- PerformaTrak

## Respironics PerformaTrak

### Key benefits

- Dual-density foam forehead pads provide patient comfort and mask stability
- With only three sizes, clinicians can quickly fit a wide range of patients
- Swivel clips simplify re-application by eliminating the need to refit after removal

Feature level



Entrainment elbow (EE)  
Also available with standard elbow (SE)

## Respironics Image3

### Key benefits

- With only three sizes, clinicians can quickly fit a wide range of patients
- Three-point headgear provides evenly distributed sealing pressure
- Swivel clips simplify re-application by eliminating the need to refit after removal

Feature level



Entrainment elbow (EE)  
Also available with standard elbow (SE)



## Respironics PerformaTrak

### Key benefits

- CapStrap headgear provides an excellent fit and simple re-application
- With only three sizes, clinicians can quickly fit a wide range of patients
- Hook-n-loop grab tabs allow quick, easy removal and adjustment of headgear

Feature level



## Less contact for more comfort

When chronically ill patients enter the hospital in respiratory distress, they typically want an NIV interface as comfortable as what they use at home. They want an unobtrusive system to maintain quality of life during their hospital stay.

### Mask selection criteria

- Oral access
- Long-term NIV

## Respironics Contour Deluxe

### Key benefits

- Dual-density foam forehead pads provide patient comfort and mask stability
- With only three sizes, clinicians can quickly fit a wide range of patients
- Hook-n-loop grab tabs allow quick, easy removal and adjustment of headgear

Feature level



## Respironics AP111

### Key benefits

- Easy sizing and fitting
- Complete kit for quick set up
- Clear field of vision

Feature level 



## Respironics PN831

### Key benefits

- Designed for pediatric patients one year or older (> 7 kg)
- Gel cushion provides a gentle seal for sensitive skin
- Softcap headgear distributes pressure evenly over a large surface area for mask stability

Feature level 



Each Respironics AP111 interface includes all four cushion sizes




## Pediatric care masks

Respironics has a complete selection of NIV masks for pediatric patients. Choose a Respironics nasal, oro-nasal, or total face mask:

- PN831 nasal  
(patients 1 year or older and > 7kg)
- AF531 oro-nasal, small  
(patients 7 years or older and > 20 kg)
- PerforMax total face, small  
(patients 7 years or older and > 20 kg)

# Patient interface selection



### Ventilator compatibility\*

Single limb circuit

Dual limb circuit

### Patient treatments\*

Bronchoscopy procedure

Aerosol medication delivery

### Patient characteristics

Claustrophobia

Facial irregularities

With no teeth

Easier mouth access

Mouth breathing

### Application indications

Rapid application

Frequent NIV patient

Post-surgical/recovery

### Headgear type

Four-point

Four-point with chin strap

CapStrap

Three-point

Bonnet

EE elbows (leak 1,2)	•	EE elbow	EE elbow
SE elbow		SE elbow	SE elbow
Bronchoscopy elbow		Bronchoscopy elbow	
		NIVO elbows	
Claustrophobia	•		
Facial irregularities	•		
With no teeth	•	•	•
Easier mouth access		with CapStrap	•
Mouth breathing	•	•	•
Rapid application	•	with CapStrap	•
Frequent NIV patient	•	•	•
Post-surgical/recovery		•	•
Four-point	•	•	
Four-point with chin strap			
CapStrap		•	•
Three-point			
Bonnet			

\* See page 22 for elbow descriptions.

EE elbow	EE elbow	•	•	•	•
SE elbow	SE elbow	•	•		
					•
					•
with CapStrap	•	•	•	•	•
	•				
with CapStrap			•		
•	•	•	•		•
•	•				
•					
					•

\*\* See page 13 for appropriate masks for pediatric patients.



# Total face masks

Respironics PerformMax multi-use		
Description	Quantity	Part number
PerformMax EE/SE, small*	1	1049037
PerformMax EE/SE, large	1	1049038
PerformMax EE/SE, extra large	1	1080423
PerformMax EE, small*	1	1047426
PerformMax EE, large	1	1047427
PerformMax EE, extra large	1	1080419
PerformMax EE, leak 2, entrainment elbow**	1	1045313
PerformMax SE standard elbow**	1	1043210
PerformMax bronchoscopy elbow**	1	1072651
PerformMax bronchoscopy elbow**	5	1072652
PerformMax small headgear	9	1045325
PerformMax large headgear	9	1045314
PerformMax extra large headgear	9	1080519
Accessory swivel replacement kit	1	7041
Nasogastric tube sealing pad	10	452040
PerformMax multi-use cleaning and disinfection guide	1	1066982



Respironics Total		
Description	Quantity	Part number
Respironics Total face mask	1	302433
Respironics Total face mask head strap	1	1003052
Respironics Total face mask head strap	9	1010609



Respironics PerformMax single-use		
Description	Quantity	Part number
PerformMax EE, small*	1	1061619
PerformMax EE, large	1	1061620
PerformMax EE, extra large	1	1080426
PerformMax EE, small*	10	1052531
PerformMax EE, large	10	1052532
PerformMax EE, extra large	5	1080424
PerformMax SE, small*	1	1061621
PerformMax SE, large	1	1061622
PerformMax SE, extra large	1	1080517
PerformMax SE, small*	10	1052533
PerformMax SE, large	10	1052534
PerformMax SE, extra large	5	1080425
PerformMax EE, leak 2, entrainment elbow**	1	1045313
PerformMax SE, standard elbow**	1	1043210
Bronchoscopy elbow**	1	1072651
Bronchoscopy elbow**	5	1072652
PerformMax small headgear	9	1045325
PerformMax large headgear	9	1045314
PerformMax extra large headgear	9	1080519
Accessory swivel replacement kit	1	7041
Nasogastric tube sealing pad	10	452040



\* Intended for use by pediatric patients 7 years or older and > 20kg

\*\* Mask not included.

## Oro-nasal masks

Respironics AF531 with CapStrap		
Description	Quantity	Part number
AF531 with EE leak 1 elbow, CapStrap headgear, small*	1	1061727
AF531 with EE leak 1 elbow, CapStrap headgear, medium	1	1061728
AF531 with EE leak 1 elbow, CapStrap headgear, large	1	1061729
AF531 with EE leak 1 elbow, CapStrap headgear, small*	10	1072627
AF531 with EE leak 1 elbow, CapStrap headgear, medium	10	1072628
AF531 with EE leak 1 elbow, CapStrap headgear, large	10	1072629
AF531 with SE standard elbow, CapStrap headgear, small*	1	1061730
AF531 with SE standard elbow, CapStrap headgear, medium	1	1061731
AF531 with SE standard elbow, CapStrap headgear, large	1	1061732
AF531 with SE standard elbow, CapStrap headgear, small*	10	1072630
AF531 with SE standard elbow, CapStrap headgear, medium	10	1072631
AF531 with SE standard elbow, CapStrap headgear, large	10	1072632
AF531 with EE/NE nebulizer elbow, CapStrap headgear, small*	10	1072618
AF531 with EE/NE nebulizer elbow, CapStrap headgear, medium	10	1072619
AF531 with EE/NE nebulizer elbow, CapStrap headgear, large	10	1072620

\* Intended for use by pediatric patients 7 years or older and > 20 kg.



Respironics AF531 with four-point headgear		
Description	Quantity	Part number
AF531 with EE leak 1 elbow, four-point headgear, small*	1	1061718
AF531 with EE leak 1 elbow, four-point headgear, medium	1	1061719
AF531 with EE leak 1 elbow, four-point headgear, large	1	1061720
AF531 with EE leak 1 elbow, four-point headgear, small*	10	1072621
AF531 with EE leak 1 elbow, four-point headgear, medium	10	1072622
AF531 with EE leak 1 elbow, four-point headgear, large	10	1072623
AF531 with EE leak 2 elbow, four-point headgear, small*	1	1061724
AF531 with EE leak 2 elbow, four-point headgear, medium	1	1061725
AF531 with EE leak 2 elbow, four-point headgear, large	1	1061726
AF531 with SE standard elbow, four-point headgear, small*	1	1061721
AF531 with SE standard elbow, four-point headgear, medium	1	1061722
AF531 with SE standard elbow, four-point headgear, large	1	1061723
AF531 with SE standard elbow, four-point headgear, small*	10	1072624
AF531 with SE standard elbow, four-point headgear, medium	10	1072625
AF531 with SE standard elbow, four-point headgear, large	10	1072626
AF531 with EE nebulizer elbow, four-point headgear, small*	10	1061733
AF531 with EE nebulizer elbow, four-point headgear, medium	10	1061734
AF531 with EE nebulizer elbow, four-point headgear, large	10	1061735

\* Intended for use by pediatric patients 7 years or older and > 20 kg.



## Oro-nasal procedure packs

Respironics AF531 CapStrap headgear procedure packs		
Description	Quantity	Part number
AF531 mask with EE leak 1 nebulizer elbow (no nebulizer), CapStrap headgear, and FEP circuit with main flow filter		
Small*	5	1084740
Medium	5	1084741
Large	5	1084742

\* Intended for use by pediatric patients 7 years or older and > 20 kg.

Note: FEP - filter exhalation port

Respironics AF531 four-point headgear procedure packs		
Description	Quantity	Part number
AF531 mask with EE leak 1 entrainment elbow, four-point headgear, and FEP circuit with main flow filter		
Small*	5	1072639
Medium	5	1072640
Large	5	1072641

\* Intended for use by pediatric patients 7 years or older and > 20 kg.

Note: FEP - filter exhalation port

# Aerogen Pro-X controllers and systems

## NIVO Pro-X controller kits\*

Description	Quantity	Part number
NIVO Pro-X controller kit – US	1	1076301
NIVO Pro-X controller kit – Canada	1	1082697

\* Each system contains: Power source, Aerogen Pro-X, equipment mount, and control module cable.

## NIVO Pro-X nebulizer systems\*

Description	Quantity	Part number
NIVO Pro-X nebulizer system – US	1	1076300
NIVO Pro-X nebulizer system – Canada	1	1082691


\* Each kit contains: Power source, Aerogen Pro-X and NIVO nebulizer; equipment mount, control module cable, and NE elbows (NE elbows have a port for NIVO nebulizer for both EE and SE elbows).

## NIVO nebulizers and elbows

Description	Quantity	Part number
NIVO nebulizer	5	1076302
NIVO nebulizer, EE, leak 1 entrainment elbow	5	1084320
NIVO nebulizer, SE standard elbow	5	1084323



# NIVO controller and system components



- |   |                       |   |
|---|-----------------------|---|
| 1 | AC/DC adapter         | Power source or charge controller                                 |
| 2 | Pro-X control module  | Can use AC/DC adapter or internal battery up to 45 minutes        |
| 3 | Equipment mount       | Use for pole or humidifier slot mounting configuration            |
| 4 | AF531 mask*           | Mask to be used   |
| 5 | Control module cable  | Connect controller to NIVO nebulizer                              |
| 6 | Nebulizer elbows (NE) | Entrainment elbow (EE) or standard elbow (SE) with nebulizer port |
| 7 | NIVO nebulizer        | Aerosol generator   |

\* NE elbow can be added to an existing AF531 or mask can be purchased with NE elbow installed.  
See page 18 for mask ordering options.

# Oro-nasal masks

## Respironics AF531 accessories and replacement parts

Description	Quantity	Part number
EE, leak 1 entrainment elbow	1	1072645
EE, leak 1 entrainment elbow	5	1072646
EE, leak 2 entrainment elbow	1	1045313
EE, leak 2 entrainment elbow	5	1072647
SE, standard elbow	1	1043210
SE, standard elbow	5	1072648
CapStrap headgear with clips	1	1072654
Four-point headgear with clips, small	1	1072655
Four-point headgear with clips, large	1	1072656
Four-point headgear with clips, extra large	5	1083188
AF531 sizing gauge	25	1078287
Swivel	1	7041
Nasogastric sealing pad	10	452040
CleanClip bracket assembly	1	1076430
CleanClip bracket thumb screw	1	1076197
Bronchoscopy elbow	1	1072651
Bronchoscopy elbow	5	1072652
EE, leak 1 nebulizing elbow (no nebulizer)	1	1084324
EE, leak 1 nebulizing elbow (no nebulizer)	5	1084325
SE, nebulizing elbow (no nebulizer)	1	1084326
SE, nebulizing elbow (no nebulizer)	5	1084327
NIVO nebulizer	5	1076302
NIVO nebulizer, EE, leak 1 entrainment elbow kit	5	1084320
NIVO nebulizer, SE standard elbow kit	5	1084323



EE, leak 1  
entrainment elbow



EE, leak 1 shown with  
NIVO nebulizer



SE, standard elbow



SE, standard elbow shown  
with NIVO nebulizer



EE, leak 2  
entrainment elbow



Bronchoscopy elbow



## Respironics AF811

Description	Quantity	Part number
AF811 oro-nasal single-use mask EE, small	1	1063047
AF811 oro-nasal single-use mask EE, medium	1	1063016
AF811 oro-nasal single-use mask EE, large	1	1063015
AF811 oro-nasal single-use mask EE, small	5	1063053
AF811 oro-nasal single-use mask EE, medium	5	1063052
AF811 oro-nasal single-use mask EE, large	5	1063051
AF811 oro-nasal single-use mask SE, small	1	1063050
AF811 oro-nasal single-use mask SE, medium	1	1063049
AF811 oro-nasal single-use mask SE, large	1	1063048
AF811 oro-nasal single-use mask SE, small	5	1063056
AF811 oro-nasal single-use mask SE, medium	5	1063055
AF811 oro-nasal single-use mask SE, large	5	1063054
CapStrap replacement headgear	5	1019543
Four-point replacement headgear, petite	5	1018472
Four-point replacement headgear, small, medium, large	5	1007365
Nasogastric tube sealing pad	10	452040

## Respironics PerformaTrak

Description	Quantity	Part number
PerformaTrak oro-nasal mask EE, small	1	1012623
PerformaTrak oro-nasal mask EE, medium	1	1012624
PerformaTrak oro-nasal mask EE, large	1	1012572
PerformaTrak oro-nasal mask EE, small	10	1018592
PerformaTrak oro-nasal mask EE, medium	10	1018593
PerformaTrak oro-nasal mask EE, large	10	1018594
PerformaTrak oro-nasal mask SE, small	1	1012636
PerformaTrak oro-nasal mask SE, medium	1	1012637
PerformaTrak oro-nasal mask SE, large	1	1012638
PerformaTrak oro-nasal mask SE, small	10	1018596
PerformaTrak oro-nasal mask SE, medium	10	1018597
PerformaTrak oro-nasal mask SE, large	10	1018598
PerformaTrak oro-nasal mask with CapStrap EE, small	1	1063833
PerformaTrak oro-nasal mask with CapStrap EE, medium	1	1063834
PerformaTrak oro-nasal mask with CapStrap EE, large	1	1063835
PerformaTrak oro-nasal mask with CapStrap SE, small	1	1063836
PerformaTrak oro-nasal mask with CapStrap SE, medium	1	1063897
PerformaTrak oro-nasal mask with CapStrap SE, large	1	1063898
PerformaTrak oro-nasal mask with CapStrap EE, small	10	1019478
PerformaTrak oro-nasal mask with CapStrap EE, medium	10	1019479
PerformaTrak oro-nasal mask with CapStrap EE, large	10	1019480
PerformaTrak oro-nasal mask with CapStrap SE, small	10	1019511
PerformaTrak oro-nasal mask with CapStrap SE, medium	10	1019525
PerformaTrak oro-nasal mask with CapStrap SE, large	10	1019533
Four-point head strap with clips	1	1015788
Four-point head strap with no clips, petite	5	1018472
Nasogastric tube sealing pad	10	452040



## Oro-nasal masks

Respironics Image3		
Description	Quantity	Part number
Image3 oro-nasal mask EE, small	1	1004878
Image3 oro-nasal mask EE, medium	1	1004849
Image3 oro-nasal mask EE, large	1	1010871
Image3 oro-nasal mask EE, small	10	1018585
Image3 oro-nasal mask EE, medium	10	1018586
Image3 oro-nasal mask EE, large	10	1018587
Image3 oro-nasal mask SE, small	1	1012526
Image3 oro-nasal mask SE, medium	1	1012528
Image3 oro-nasal mask SE, large	1	1012529
Image3 oro-nasal mask SE, small	10	1018578
Image3 oro-nasal mask SE, medium	10	1018579
Image3 oro-nasal mask SE, large	10	1018580
Image3 standard replacement headgear	1	1004877
Image3 Softcap replacement headgear	1	1006416



## Nasal masks

Respironics PerformaTrak		
Description	Quantity	Part number
PerformaTrak nasal mask, CapStrap headgear, petite	10	1048709
PerformaTrak nasal mask, CapStrap headgear, small	10	1048708
PerformaTrak nasal mask, CapStrap headgear, medium/ large	10	1048710
Nasogastric tube sealing pad	10	452040



Respironics Contour Deluxe		
Description	Quantity	Part number
Contour Deluxe nasal mask, petite	5	1016980
Contour Deluxe nasal mask, small	5	1016693
Contour Deluxe nasal mask, medium/large	5	1016691
Contour Deluxe nasal mask, petite	10	1016982
Contour Deluxe nasal mask, small	10	1016694
Contour Deluxe nasal mask, medium/large	10	1016692
Contour Deluxe head strap replacement, petite	5	1018472
Contour Deluxe head strap replacement, small and medium/large	5	1007365



## Nasal masks

### Respironics AP111

Description	Quantity	Part number
AP111 interface with cradles*	1	1062128
AP111 interface with cradles*, circuit, and main flow filter	1	1062130
AP111 interface with cradles*, circuit, and flow filter	5	1062132
AP111 headgear	5	1062134

\*All four cushion sizes included with each interface frame.

### Respironics PN831\*

Description	Quantity	Part number
PN831 with Softcap	1	1072201
PN831 with Softcap, circuit, and main flow filter	1	1072202
PN831 with Softcap	5	1072203
PN831 with Softcap, circuit and main flow filter	5	1072204
Softcap replacement headgear, child	1	302465

\* Intended for use by pediatric patients 1 year or older and > 7 kg.



## Procedure packs

### Procedure packs

Description	Quantity	Part number
AF531 mask with EE leak 1 nebulizer elbow (no nebulizer), CapStrap headgear, and FEP circuit with main flow filter		
Small*	5	1084740
Medium	5	1084741
Large	5	1084742
AF531 mask with EE leak 1 entrainment elbow, four-point headgear, and FEP circuit with main flow filter		
Small*	5	1072639
Medium	5	1072640
Large	5	1072641
AF811 gel mask with EE leak 1 entrainment elbow, CapStrap headgear, and FEP circuit with main flow filter and exhalation filter		
Small	5	1074839
Medium	5	1074838
Large	5	1074837
AF811 gel mask with EE leak 1 entrainment elbow, CapStrap headgear, and FEP circuit with main flow filter		
Small	5	1074061
Medium	5	1074060
Large	5	1074059
PerforaTrak mask with EE leak 1 entrainment elbow, four-point headgear, and FEP circuit with main flow filter		
Small	5	1074055
Medium	5	1074054
Large	5	1074053
PerforaTrak mask with EE leak 1 entrainment elbow, CapStrap headgear, and FEP circuit with main flow filter		
Small	5	1074058
Medium	5	1074057
Large	5	1074056
PerforMax (single use) mask with EE leak 1 entrainment elbow and FEP circuit with main flow filter		
Small	5	1074063
Large	5	1074062

\* Intended for use by pediatric patients 7 years or older and > 20 kg.

Note: FEP - filter exhalation port


# Circuits and connections

Respironics supports its noninvasive ventilation systems with a variety of single-use circuits, including heated wire and non-heated-wire.



Circuit	Description	Part number	Quantity	V60	BiPAP Vision	Trilogy 202	BiPAP Focus	BiPAP S/T-D
1	Invasive circuit with exhalation port, water trap and temperature probe ports, proximal airway filter and humidifier coupling tube, hanger, and hose clips	652002	10	•	•	•		
1	Invasive circuit with exhalation port, water trap and temperature probe ports, proximal airway filter and humidifier coupling tube, hanger, and hose clips	652001	20	•	•	•		
2	Single limb (25mm) heated wire circuit with exhalation port	1020523	10	•	•	•		
3	Noninvasive circuit with main flow filter, disposable swivel exhalation port with filter port and filter, proximal pressure line, hanger, and hose clips	1065830	10	•	•	•	•	
4	Noninvasive circuit with main flow filter, disposable swivel exhalation port with filter port, proximal pressure line, hanger, and hose clips (exhalation filter not included)	1065832	10	•	•	•	•	
5	Noninvasive circuit with disposable swivel exhalation port with filter port, proximal pressure line, hanger, and hose clips (main flow and exhalation filters not included)	1069210	10	•	•	•	•	
6	Noninvasive circuit with disposable exhalation port, proximal pressure line, hanger, and hose clips	582073	10	•	•	•	•	
7	BiPAP S/T-D circuit with exhalation port, proximal pressure line, hanger, and hose clips	312107	10					•
*	Trilogy adult invasive passive disposable non-heated wire circuit	1069443	10			•		
*	Trilogy adult invasive passive disposable non-heated wire circuit	1069444	20			•		
*	Trilogy adult noninvasive passive disposable non-heated wire circuit	1069467	10			•	•	
*	Trilogy adult noninvasive passive disposable non-heated wire circuit	1069468	20			•	•	




\* Circuit not shown



# Mask components

Illustration	Components	Quantity	Part number
1	Plateau exhalation valve	1	302312
2	Plateau exhalation valve diaphragm	5	302310
3	Whisper Swivel II	1	332113
4	Disposable exhalation port with filter exhalation port (no filter)	10	1065775
4a	Replacement port caps	5	1061825
5	Disposable exhalation port	10	312149
6	Oxygen enrichment attachment	1	312710
6	Oxygen enrichment attachment	10	312010
7	Bacteria/virus filter	1	1014047
7	Bacteria/virus filter	10	342077
8	Proximal airway filter	1	1002362
9	Tube management clip	2	1039967
10	EE, leak 1 entrainment elbow replacement kit	1	1072645
10	EE, leak 1 entrainment elbow replacement kit	5	1072646
11	EE, leak 2 entrainment elbow replacement kit	1	1045313
11	EE, leak 2 entrainment elbow replacement kit	5	1072647
12	SE, standard elbow	1	1043210
12	SE, standard elbow	5	1072648
13	Bronchoscopy elbow	1	1072651
13	Bronchoscopy elbow	5	1072652
14	Swivel replacement kit (22 mm male-male)	1	7041
15	CleanClip bracket assembly	1	1076430
15	CleanClip bracket thumb screw	1	1076197
16	Nasogastric tube sealing pad	10	452040
17	PerforMax sizing gauge, small, large and extra large	1	1095097
18	Replacement ball and socket	2	1009723
	PerformaTrak sizing gauge (not shown)	1	1047937
*	EE, leak 1 nebulizing elbow (no nebulizer)	1	1084324
*	EE, leak 1 nebulizing elbow (no nebulizer)	5	1084325
*	SE nebulizing elbow (no nebulizer)	1	1084326
*	SE nebulizing elbow (no nebulizer)	5	1084327
*	NIVO nebulizer	5	1076302
*	NIVO nebulizer, EE, leak 1 entrainment elbow	5	1084320
*	NIVO nebulizer, SE standard elbow	5	1084323

\* See page 21 for illustrations.







# Headgear

Illustration	Headgear	Size	Compatibility	Quantity	Part number
1	CapStrap headgear with ball and socket clips	Adult (one size)	AF811, PerformaTrak	5	1019543
1	CapStrap headgear with clips	Adult (one size)	AF531	1	1072654
2	Three-point Softcap headgear	Adult (one size)	Image3	1	1006416
3	Three-point headgear	Adult (one size)	Image3	1	1004877
4	Four-point headgear with ball and socket clips	Adult (one size)	AF811, PerformaTrak	1	1015788
5	Four-point headgear, no clips	Petite	AF811, PerformaTrak	5	1018472
			Contour Deluxe		
5	Four-point headgear, no clips	Adult (one size)	AF811, PerformaTrak	5	1007365
			Contour Deluxe		
6	Four-point headgear	Adult (one size)	Total	1	1003052
6	Four-point headgear	Adult (one size)	Total	9	1010609
7	Four-point headgear	Adult (one size)	Vinyl nasal	10	312106
8	Four-point headgear with crown strap, clips	Small	PerforMax	5	1051857
8	Four-point headgear with crown strap, clips	Large	PerforMax	5	1051816
8	Four-point headgear with crown strap, clips	Extra large	PerforMax	5	1080518
8	Four-point headgear with crown strap, clips	Small	PerforMax	9	1045325
8	Four-point headgear with crown strap, clips	Large	PerforMax	9	1045314
8	Four-point headgear with crown strap, clips	Extra large	PerforMax	9	1080519
8	Four-point headgear with crown strap, clips	Small	AF531	1	1072655
8	Four-point headgear with crown strap, clips	Large	AF531	1	1072656
8	Four-point headgear with crown strap, clips	Extra large	AF531	5	1083188
9	Two-point headgear with chin strap	Adult (one size)	AP111	5	1062134


1




2




3




4




5




6




7



8



9



**Philips Healthcare is part of  
Royal Philips Electronics**

How to reach us

[www.philips.com/healthcare](http://www.philips.com/healthcare)

[healthcare@philips.com](mailto:healthcare@philips.com)

Asia

+49 7031 463 2254

Europe, Middle East, Africa

+49 7031 463 2254

Latin America

+55 11 2125 0744

North America

+1 425 487 7000

800 285 5585 (toll free, US only)

Please visit [www.philips.com/NIVmasks](http://www.philips.com/NIVmasks)



© 2012 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands.  
4522 962 80371 \* JAN 2012