



# iQuest

STANDARD • HOSPITAL • LOUNGE





In 10 years, the Quest examination couch has satisfied more than 15,000 healthcare professionals. General practitioners, cardiologists, dermatologists and physiotherapists from over 70 countries have come to appreciate

the French quality of this robust, comfortable and easy to access examination couch. Five years ago, Promotal conducted the largest survey ever performed in Europe and the Middle East. This 3-year study, carried out in collaboration with more than 200 health professionals, revealed 3 basic needs: equal treatment of people with mobility limitations, efficiency in carrying out medical examinations and user-friendly equipment. iQuest is Promotal's new revolution. We have designed an examining couch with a sleek, lightweight design that will improve your comfort and that of your patients thanks to the low position that is the lowest on the market, the programmable positions and the autonomous movements that are 2 times faster.



## EQUAL TREATMENT

### MAKE THINGS EASIER FOR YOUR PATIENTS WITH MOBILITY LIMITATIONS :

- > Low 46 cm position provides easier access for patients
- > Patients can use the raised backrest for support when getting on the couch
- > Safe, smooth electrical movements

### EXAMINE ALL YOUR PATIENTS EASILY :

- > Accommodates up to 250 Kg
- > Superior comfort with the Lounge upholstery



## EFFICIENT DESIGN

### YOUR EXAMINATIONS WILL GO SMOOTHLY THANKS TO

- > Movements that are 2 times faster
- > The intuitive position of the patient on the Lounge upholstery
- > The autonomous movements: Press once to change positions
- > The intuitive controls and the pre-programmed patient access position



## USER-FRIENDLY DESIGN

### CHANGE POSITIONS EFFORTLESSLY :

- > All-electric examination couch: Height, backrest, Trendelenburg / reverse Trendelenburg
- > Variable height from 46 to 95 cm
- > 2 programmable positions, 100% configurable
- > Intuitive controls built into the backrest (Lounge model)

### IMPROVE YOUR PATIENTS' COMFORT THANKS TO :

- > The superior comfort Lounge upholstery



ACCESS



EXAMINATION



HIGH  
POSITION



TRENDLENBURG



REVERSE  
TRENDLENBURG



3 MODELS,  
3 NEEDS.

● OPTIONS  
● ACCESSORIES



**iQuest**  
STANDARD

Ref. 3050-40  
With retractable wheels  
Ref. 3050-50  
Upholstery width 65 cm  
Ref. 600150-65  
Upholstery width 75cm  
Ref. 600150-75



**iQuest**  
LOUNGE  
Ref. 3055-20



**iQuest**  
HOSPITAL  
Ref. 3050-60




Face cut-out  
65cm Ref. 2231-65  
75cm Ref. 2231-75




Factory-installed  
upholstery  
Ref. 600000-01



Battery  
Ref. 3092-01



Circuit breaker  
Ref. 3058-01



Cable reel  
Ref. 3032-01XX



Anti-static wheels  
Ref. 3031-01




Integrated plug  
EU-US-UK Ref. 30400-01



Equipotential plug  
Ref. 2059-01



Clamp  
Ref. 166



Pair of clamps on rail  
Ref. 878-12



Pair of legrests  
Ref. 840



Pair of legrest covers  
Ref. 861



Pair of stirrups  
Ref. 850



IV pole with 2 hooks  
Ref. 985-01




Auto-locking IV pole  
with 2 hooks  
Ref. 2985-01




Armrest  
(one)  
Ref. 2536



Armrest cover  
(one)  
Ref. 2537




Lounge headrest  
Ref. 30679-01



Flat upholstery headrest  
Ref. 1102-10



Adjustable paper roll  
holder  
Ref. 2056-02



Backrest extension  
Ref. 2055-01



Wired foot control  
Ref. 3051-10



Plastic cover (foot end)  
for flat upholstery  
65cm Ref. 600150-65C  
75cm Ref. 600150-75C



Plastic cover  
(head end)  
Ref. 3056-02



Plastic cover  
(foot end)  
Ref. 3056-01




Pair of barriers  
Ref. 85114-01




Pair of head rails  
Ref. 30878-XXT




Pair of foot rails  
Ref. 30878-XXP



STANDARD



HOSPITAL



LOUNGE

32 - BLUE		CLASSIC
93 - COTTAGE		
35 - RED		
34 - YELLOW		
91 - WHITE		
92 - MARINE		VITAMIN
56 - DARKNESS		
57 - ANTISTATIC		
64 - GRAPE		
67 - LAGOON		
66 - SWEET		COSY
63 - CEDRAT		
80 - SKY		
81 - SEA GREEN		
65 - LILAC		
96 - TIN		
82 - MOUSE GREY		
83 - NATURAL		
84 - CHOCOLATE		