



SOLID. RELIABLE. SAFE.

# Cart light duty

Mobile solutions for every application



## CART LIGHT DUTY

A roll stand for light applications of up to 14 kg/ 30.9 lbs. All carts are tested and certified as safe before supply. The robust construction makes it a favourite amount nursing staff.

### DESIGNED FOR


- + Daily repetitive use
- + Reliability
- + Current regulation (EU) 2017/745 (MDR) & IEC 60601-1 safety standards
- + FDA registered and 510k compliant

### FEATURES

- + Easy to clean
- + No pooling areas for liquids
- + Easy-glide movement
- + Lockable castors for stability
- + Small footprint
- + 5-year warranty


### ACCEPTS

- + All standard fixings
- + External power supply holsters
- + Baskets
- + Many other accessories




## MONITOR PRESENTATIONS

Which monitor would you like to mount?  
 CIM med offers solutions for a wide range of medical device manufacturers.



## ACCESSORIES

A huge variety of accessories is available which you will find in our accessories catalogue.  
 Please note that not every accessory is suitable for every cart.



# Cart light duty - for 5" plunge plate applications

A



Illustration shown with Dräger Infinity Delta

<b>Monitor adaption</b>	5" plunge plate with tilt function
<b>Accessories</b>	5-star aluminum base with counter weight double wheel castors, ø 100 mm / 3.9", 2 lockable column to accept components and accessories on opposite sides
<b>Max. load capacity</b>	14 kg / 30.9 lbs
<b>Optional accessories</b>	see page 3


### Suitable for

- Generic 5" plunge plate solutions
- Dräger monitoring: Infinity Delta
- Fukuda Denshi monitoring: Dynascope DS-8007
- GE monitoring: Procure B20/ B40, GE Carescape B105/ B125/ B450/ B650, V 100/ VC 100
- Mindray monitoring: BeneVision N15/ N17, BeneVision N12/ iMEC/ uMEC/ iPM/ ePM, BeneView T1/ N1/ TDS/ T5/ T8, T9, VS-600, VS-900
- Nihon Kohden monitoring: Life Scope PT BSM 1700/ CMS 1500/ TR BSM 6000, VS BSM 3500/ 3700, Vismo PVM 2700
- Philips monitoring: Efficia CM100/ CM120/ CM150, Sure Signs/ Vital Signs VS2+/ VSi/ VS4

Please note that monitor-specific mounting adapter are required for certain manufacturers.

## Cart light duty - for VESA 75/100 applications

K



### Monitor adaption

CL-bracket, VESA 75/100

### Accessories

5-star aluminum base with counter weight  
double wheel castors,  $\varnothing$  100 mm / 3.9", 2 lockable  
column to accept components and accessories on opposite sides

### Max. load capacity

10 kg / 22 lbs

### Optional accessories

see page 3


### Suitable for

- Computer hardware with VESA 75/100
- Dräger monitoring: Infinity C500/700, Explorer
- GE monitoring: Carescape B850
- Nihon Kohden monitoring: Life Scope G7/ G9
- Spacelabs monitoring: Xprezzon display

All rear mounted flat panels of various manufacturers.

## Cart light duty - for notebooks

TN



**Monitor adaption**

360 x 280 mm / 14.2 x 11", with tilt function

**Accessories**

5-star aluminum base without counter weight  
 double wheel castors,  $\varnothing$  100 mm / 3.9", 2 lockable  
 column to accept components and accessories on opposite sides

**Max. load capacity**

7 kg / 15.4 lbs (without counter weight)

**Optional accessories**

see page 3



**CIM med GmbH**

Euro-Industriepark  
Margot-Kalinke-Str. 9  
80939 Munich / Germany

T +49 89. 978 94 08-0

---

**CIM med Inc.**

2024 E Independence Commerce Dr  
Matthews, NC 28105 / USA

T +1 855 895 4722

---

[info@cim-med.com](mailto:info@cim-med.com)  
[www.cim-med.com](http://www.cim-med.com)

