



Rotator “Orbit-II”

FEATURES

- Adjustable speed from 10 to 80 r.p.m.
- Adjustable device for the inclination of the rotation axis.
- 9 configurable programmes.
- IP 21 protection.
- Brushless induction motor.
- Quick operation button of 30 seconds.
- Programmable time from 1 to 99 minutes.
- Recommended for cell and haematology extractions and cultures.
- Epoxy coated outer casing.

MODEL

Code	Speed	Height / Width / Depth			Consumption	Weight
3000446	10 a 80	38	32	32	4.5	4.5



ACCESSORIES

Adapters for Erlenmeyers, with stainless steel clamps

Code	Capacity
3000850	8 Erlenmeyers of 50 ml.
3006100	6 Erlenmeyers of 100 ml.
3005250	5 Erlenmeyers of 250 ml.

Tube adapters, aluminium, with P.V.C. coated clamps for different tube diameters.

Code	Capacity	Code	Capacity
3000356	35 tubes of 6 to 9 mm Ø.	3002215	22 tubes of 17 to 20 mm Ø.
3003011	30 tubes of 10 to 13 mm Ø.	3001820	18 tubes of 21 to 24 mm Ø.
3002512	25 tubes of 14 to 16 mm Ø.	3001723	17 tubes of 23 to 28 mm Ø.



Tilting Mixer Shaker “Movil-Tub”

FEATURES

Fixed speed tilting shaker mixer with removable silicon, anti slide tray. Recommended especially for haematology, chemistry and microbiology applications.

MODEL

Part No.	Tube capacity	Tube Ø / length mm	Speed r.p.m.	Height / Width / Depth (exterior) cm			Tilt angle	Power W	Weight Kg
7000913	12	20 160	15	12	34	20	25°	6	2,5



Roller and tilt mixer “Movil-Rod”

WITH VARIABLE SPEED CONTROL FROM 10 TO 80 r.p.m.

FEATURES

Roll and tilt mixer with 7 rollers, double action, gently rolls the tubes while rocking, produces a sinusoidal movement, providing excellent sample homogenisation. Recommended for suspensions of biological and haematological fluids.

MODEL

Part No.	For tube Ø up to mm	Speed r.p.m.	Length (rollers) cm	Height / Width / Depth (exterior) cm			Power W	Weight Kg
7001723	35	10 - 80	24	9.8	34.9	26.5	12	3.6





Roller and tilt mixer “Movil-Rod II”

FIXED SPEED 36 RPM.

FEATURES

Roll and tilt mixer with 5 rollers, double action, gently rolls the tubes while rocking, produces a sinusoidal movement, providing excellent sample homogenisation. Recommended for suspensions of biological and haematological fluids.

MODEL

Part No.	For tube Ø up to mm	Speed r.p.m.	Length (rollers) cm	Height / Width / Depth (exterior) cm			Power W	Weight Kg
7001726	35	36	23	14	26	15,5	6	3



Rocking Mixer “Vibromatic”

DIGITAL ELECTRONIC CONTROL OF SPEED AND TIME.



Optimum shaking frequency: 950 oscillations per minute

APPLICATIONS

Dissolution, extractions and homogenisation.

FEATURES

Electronic control from 100 to 950 oscillations per minute. (1900 vibrations) and timer from 1 to 999 minutes or continuous operation.

Oscillation amplitude: 8 mm.

Easily detachable arms equipped with 8 clamps that allow all types of tubes and flask with diameters of between 5 to 50 mm can be held in position.

Maximum load: 4 kg.


CONTROL PANEL

1. Digital display indicator of osc./min. and time.
2. Parameter display selector.
3. Parameter selector, speed or time.
4. Push button increase value.
5. Push button decrease value.
6. Push button start-stop.



MODEL

Part No.	Oscillations p.m.	Maximum capacity Kg	Height / Width / Depth (exterior) cm			Power W	Weight Kg
7000384	100 - 950	4	23	77	27	34	11



ACCESSORIES

Platform support for Erlenmeyer flasks, made from AISI 304 stainless steel with clip adapter position holes for 50 ml to 1000 ml flasks. Platform dimensions: 38 cm long x 19 cm wide.

Weight: 2 Kg.

Platform capacity:

- 18 Erlenmeyer flasks of 50 ml.
- 12 Erlenmeyer flasks of 100 ml.
- 8 Erlenmeyer flasks of 250 ml.
- 6 Erlenmeyer flasks of 500 ml.
- 3 Erlenmeyer flasks of 1000 ml.

Part No. **3001221**



Model Vibromatic with adapter platform.

Adapters for Erlenmeyer flasks.

Made from tempered stainless steel.

Part No.

- | | | |
|---------|-------------|----------|
| 1001003 | Adapter for | 50 ml. |
| 1000978 | Adapter for | 100 ml. |
| 1000979 | Adapter for | 250 ml. |
| 1000980 | Adapter for | 500 ml. |
| 1000981 | Adapter for | 1000 ml. |

