



BioFusion P Series

P300 Syringe Pump

P300 / P300 Neo / P300 Ex / P300 Pro

Guangdong Biolight Meditech Co.,Ltd.

Add: NO.2 Innovation First Road, Technical Innovation Coast, Hi-tech Zone, Zhuhai, P.R.China,
Postcode: 519085 Tel: +86-756-3399935 Fax: +86-756-3399911 E-Mail: overseas@blt.com.cn Web: global.blt.com.cn

BioFusion P300 series V1 - Aug 2025

*Specifications subject to change without prior notice.
P 300 CY - 202512



Safety



Efficient



Network



Drug library

Precision for Advanced Therapies

Designed to deliver precision, safety, and versatility, the P300 provides accurate drug infusion for applications ranging from routine care to advanced anesthesia management.

Its intuitive operation, advanced therapy modes, and compact form factor make it an ideal choice for both bedside and mobile clinical settings.



Targeted Accuracy



Clinical Safety



Comprehensive Modes



Workflow Efficiency


Model	BioFusion P300 /P300 Neo /P300 Ex /P300 Pro
Dimension	265mm x 160mm x 75mm (LxWxH)
Display	6.9 inches capacitive touchscreen
Resolution	280 x 1424 pixels

Streamlined for Every Clinician

- Easy to use, precise in delivery.
- Simple setup with guided touchscreen.
- Clear menus focus on patient care.



Enhanced Safety & Intelligent Delivery



smart infusion

- Ensures precise safety monitoring.
- Reduces unnecessary care interruptions.
- Provides prompt, customizable alerts.

Target-Controlled Infusion (TCI)


- Personalized drug delivery for anesthesia.
- Uses advanced models for stable dosing.
- Ensures optimal patient outcomes.

When managing post-operative pain, the **Patient-Controlled Analgesia (PCA)** mode empowers patients to self-administer bolus doses within clinician-defined limits. PCA offers both safety and comfort, improving the overall patient experience.

Drug Management for Injection

Drug Library

- Includes built-in safety limits.
- Offers quick preset protocols.
- Supports specialized medication delivery.



Fluid Balance Monitoring

BioFluid Solution

- Connects devices for live tracking.
- Accurately monitors fluid balance.
- Supports critical care decisions.

