

FOCUSED ON C-ARM

SPECIAL-ORTHO

Built-in support for extremity
orthopedic surgery



SPECIAL-ORTHO

Built-in support for extremity ortopedic surgery

This built-in support is a Simad Moonray System accessory designed and realized to make the C-Arm effective for ortopedic application on extremities.

It replaces the standard support, usually provided of telescopic support and connected to the table.

It can be installed on the I.I. and assure a better positioning of the limb in terms of image quality, field of view and x-ray scattering reduction.


Being the Special-Ortho support adjustable, and being it without any structure to be hooked at the floor or table, it allows the user to come near the patient in the easiest way without any physical impedance.

Also, this extremities support increases the safety level of the surgical activity since it represents an additional protection to avoid accidental damages.

This support is easily installable and removable and fits with all Moonray models.

It is made of a particular material which makes it particularly radio-transparent, in order to prevent from any increase of the x-ray dose to get an adequate imaging.

For particular kind of procedures, the Special-Ortho support can be provided of hand traction accessory.



X-ray scattering cutback

Largest dimensions and superior image quality

Perfect for extremity surgical applicatons



www.simad.net

SIMAD s.r.l.
Società Unipersonale

"Direzione e coordinamento di
Sago Medica s.r.l."

Via XXV Aprile, 22
41037 Mirandola - Modena - Italy
Tel. +39 0535 24 305
Fax: +39 0535 25 010
info@simad.net

