

The Lapro-Flex[®] Family of Articulating Retractors




Mediflex's family of Lapro-Flex[®] articulating laparoscopic retractors allows for convenient and atraumatic retraction of organs and tissue in the peritoneal space.

Mediflex[®]
SURGICAL PRODUCTS

LF-1113

The Lapro-Flex® Self-forming Articulating Retractor

These retractors can be introduced through a 5mm cannula and articulated to their shape by rotating a knob in the handle. They are available in several different configurations including right angle, triangular and angled triangular.




Ref. No.	Description	Shaft Length
91680	90° Articulation, Finger-Styled	38cm
91682	60mm Triangle Base	32cm
91682-A*	60mm Angled Triangle Base	32cm
91683	80mm Triangle Base	32cm
91683-A*	80mm Angled Triangle Base	32cm



Ref. No. 91680

Ref. No. 91682
&
Ref. No. 91682-A

Ref. No. 91683
&
Ref. No. 91683-A



The rigid arm shown on the other side is Mediflex product 73000-DT

For Additional Information Contact Mediflex
or your Local Mediflex Distributor



Mediflex®
SURGICAL PRODUCTS

www.mediflex.com
Toll-free: 1-800-764-0474
Phone: 1-631-582-6424
Fax: 1-631-582-4520

© 2013 All rights reserved. Lapro-Flex and Mediflex are registered trademarks of Flexbar Machine Corporation.